



PGAM5 Antibody, Biotin conjugated

Product Code	CSB-PA839338LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HS1
Immunogen	Recombinant Human Serine/threonine-protein phosphatase PGAM5, mitochondrial protein (30-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein phosphatase PGAM5, mitochondrial (EC 3.1.3.16) (Bcl-XL-binding protein v68) (Phosphoglycerate mutase family member 5), PGAM5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PGAM5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.