



NPL Antibody, Biotin conjugated

Product Code	CSB-PA874830LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXD5
Immunogen	Recombinant Human N-acetylneuraminate lyase protein (17-77AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-acetylneuraminate lyase (NALase) (EC 4.1.3.3) (N-acetylneuraminate pyruvate-lyase) (N-acetylneuraminic acid aldolase) (Sialate lyase) (Sialate-pyruvate lyase) (Sialic acid aldolase) (Sialic acid lyase), NPL, C1orf13
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NPL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.