



NPL Antibody, Biotin conjugated

Product Code	CSB-PA874830LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXD5
Immunogen	Recombinant Human N-acetylneuraminate lyase protein (17-77AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-acetylneuraminate lyase (NALase) (EC 4.1.3.3) (N-acetylneuraminate pyruvate-lyase) (N-acetylneuraminic acid aldolase) (Sialate lyase) (Sialate-pyruvate lyase) (Sialic acid aldolase) (Sialic acid lyase), NPL, C1orf13
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NPL