



BANF1 Antibody, FITC conjugated

Product Code	CSB-PA002550LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75531
Immunogen	Recombinant Human Barrier-to-autointegration factor protein (2-89AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Barrier-to-autointegration factor (Breakpoint cluster region protein 1) [Cleaved into: Barrier-to-autointegration factor, N-terminally processed], BANF1, BAF BCRG1
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	BANF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.