



CD79A Antibody, HRP conjugated

Product Code	CSB-PA004957LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11912
Immunogen	Recombinant Human B-cell antigen receptor complex-associated protein alpha chain protein (70-142AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	B-cell antigen receptor complex-associated protein alpha chain (Ig-alpha) (MB-1 membrane glycoprotein) (Membrane-bound immunoglobulin-associated protein) (Surface IgM-associated protein) (CD antigen CD79a), CD79A, IGA MB1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD79A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.