



CD79A Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA004957LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P11912 |
| Immunogen | Recombinant Human B-cell antigen receptor complex-associated protein alpha chain protein (70-142AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | B-cell antigen receptor complex-associated protein alpha chain (Ig-alpha) (MB-1 membrane glycoprotein) (Membrane-bound immunoglobulin-associated protein) (Surface IgM-associated protein) (CD antigen CD79a), CD79A, IGA MB1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | CD79A |