



CNN3 Antibody, Biotin conjugated

Product Code	CSB-PA005657LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15417
Immunogen	Recombinant Human Calponin-3 protein (2-329AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calponin-3 (Calponin, acidic isoform), CNN3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CNN3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.