



isl1 Antibody, HRP conjugated

Product Code	CSB-PA011846LB01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53405
Immunogen	Recombinant Zebrafish Insulin gene enhancer protein isl-1 protein (1-349AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Insulin gene enhancer protein isl-1 (Islet-1), isl1, isl-1
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Neuroscience
Target Names	isl1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.