



KIF11 Antibody, HRP conjugated

Product Code	CSB-PA012307LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52732
Immunogen	Recombinant Human Kinesin-like protein KIF11 protein (920-1038AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF11 (Kinesin-like protein 1) (Kinesin-like spindle protein HKSP) (Kinesin-related motor protein Eg5) (Thyroid receptor-interacting protein 5) (TR-interacting protein 5) (TRIP-5), KIF11, EG5 KNSL1 TRIP5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KIF11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.