



P4HA1 Antibody, FITC conjugated

Product Code	CSB-PA017339LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13674
Immunogen	Recombinant Human Prolyl 4-hydroxylase subunit alpha-1 protein (256-525AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prolyl 4-hydroxylase subunit alpha-1 (4-PH alpha-1) (EC 1.14.11.2) (Procollagen-proline,2-oxoglutarate-4-dioxygenase subunit alpha-1), P4HA1, P4HA
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	P4HA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.