



PRKCZ Antibody, FITC conjugated

Product Code	CSB-PA15889C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05513
Immunogen	Recombinant Human Protein kinase C zeta type protein (494-575AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein kinase C zeta type (EC 2.7.11.13) (nPKC-zeta), PRKCZ, PKC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRKCZ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.