



Sporamin A Antibody, HRP conjugated

Product Code	CSB-PA320416HB01INB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14715
Immunogen	Recombinant Ipomoea batatas Sporamin A protein (24-219AA)
Raised In	Rabbit
Species Reactivity	Ipomoea batatas
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sporamin A (Clone PIM0335)
Immunogen Species	Ipomoea batatas (Sweet potato) (Convolvulus batatas)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.