



csrA Antibody, HRP conjugated

Product Code	CSB-PA543019HB01ENT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B1XCM4
Immunogen	Recombinant Escherichia coli Translational regulator CsrA protein (1-61AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Translational regulator CsrA (Carbon storage regulator), csrA
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	csrA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.