




GRIK4 Antibody

Product Code	CSB-PA619054LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16099
Immunogen	Recombinant Human Glutamate receptor ionotropic, kainate 4 protein (826-956AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate receptor ionotropic, kainate 4 (GluK4) (Excitatory amino acid receptor 1) (EAA1) (Glutamate receptor KA-1) (KA1), GRIK4, GRIK
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRIK4

Image



Western Blot

Positive WB detected in: MCF-7 whole cell lysate

All lanes: GRIK4 antibody at 2.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 108 kDa

Observed band size: 97 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.