



FHL5 Antibody, FITC conjugated

Product Code	CSB-PA725899LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TD97
Immunogen	Recombinant Human Four and a half LIM domains protein 5 protein (16-284AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Four and a half LIM domains protein 5 (FHL-5) (Activator of cAMP-responsive element modulator in testis) (Activator of CREM in testis), FHL5, ACT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FHL5