



DNAH5 Antibody, Biotin conjugated

Product Code	CSB-PA819476LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TE73
Immunogen	Recombinant Human Dynein heavy chain 5, axonemal protein (2364-2560AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dynein heavy chain 5, axonemal (Axonemal beta dynein heavy chain 5) (Ciliary dynein heavy chain 5), DNAH5, DNAHC5 HL1 KIAA1603
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DNAH5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.