



VPS33A Antibody, Biotin conjugated

Product Code	CSB-PA850260LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AX1
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 33A protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 33A (hVPS33A), VPS33A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VPS33A