



YIPF3 Antibody, Biotin conjugated

Product Code	CSB-PA863928LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZM5
Immunogen	Recombinant Human Protein YIPF3 protein (2-148AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein YIPF3 (Killer lineage protein 1) (Natural killer cell-specific antigen KLIP1) (YIP1 family member 3) [Cleaved into: Protein YIPF3, 36 kDa form III], YIPF3, C6orf109 KLIP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	YIPF3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.