



TIMP1 Antibody, Biotin conjugated

Product Code	CSB-PA023560YD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01033
Immunogen	Recombinant Human Metalloproteinase inhibitor 1 protein (24-207AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metalloproteinase inhibitor 1 (Erythroid-potentiating activity) (EPA) (Fibroblast collagenase inhibitor) (Collagenase inhibitor) (Tissue inhibitor of metalloproteinases 1) (TIMP-1), TIMP1, CLGI TIMP
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TIMP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.