



FTL Antibody Pair

Product Code	CSB-EAP114091
Species Reactivity	Human
Tested Applications	S-ELISA
Form	Liquid
Alias	Ferritin light chain,Ferritin L subunit,FTL
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	FTL
Protein Names	Ferritin light chain
Notes	We recommend using the capture antibody at a concentration of 1ug/ml and the detection antibody at a concentration of 1ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Mouse Detection: Mouse
Components	Capture: CSB-EAP114091C Detection: CSB-EAP114091D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4