




ETV7 Antibody

Product Code	CSB-PA897576LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y603
Immunogen	Recombinant Human Transcription factor ETV7 protein (134-218AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor ETV7 (ETS translocation variant 7) (ETS-related protein Tel2) (Tel-related Ets factor) (Transcription factor Tel-2), ETV7, TEL2 TELB TREF
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ETV7

Image



Western Blot

Positive WB detected in: Rat liver tissue


All lanes: ETV7 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 39, 33, 31, 37, 30, 23 kDa

Observed band size: 39 kDa



Immunohistochemistry of paraffin-embedded human cervical cancer using CSB-PA897576LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.