



A2ML1 Antibody, Biotin conjugated

Product Code	CSB-PA001005LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8K2U0
Immunogen	Recombinant Human Alpha-2-macroglobulin-like protein 1 protein (625-751AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-2-macroglobulin-like protein 1 (C3 and PZP-like alpha-2-macroglobulin domain-containing protein 9), A2ML1, CPAMD9
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	A2ML1