



ACAT1 Antibody, Biotin conjugated

Product Code	CSB-PA001134LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24752
Immunogen	Recombinant Human Acetyl-CoA acetyltransferase, mitochondrial protein (76-106AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acetyl-CoA acetyltransferase, mitochondrial (EC 2.3.1.9) (Acetoacetyl-CoA thiolase) (T2), ACAT1, ACAT MAT
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ACAT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.