





ACTR1A Antibody

Product Code	CSB-PA001246EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61163
Immunogen	Recombinant Human Alpha-centractin protein (187-376AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-centractin (Centractin) (ARP1) (Actin-RPV) (Centrosome-associated actin homolog), ACTR1A, CTRN1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ACTR1A

Image



Western Blot
 Positive WB detected in: MCF-7 whole cell lysate
 All lanes: ACTR1A antibody at 1µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 43 kDa
 Observed band size: 43 kDa



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA001246EA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.