



ACTR1A Antibody, HRP conjugated

Product Code	CSB-PA001246EB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61163
Immunogen	Recombinant Human Alpha-centractin protein (187-376AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-centractin (Centractin) (ARP1) (Actin-RPV) (Centrosome-associated actin homolog), ACTR1A, CTRN1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ACTR1A