



ALDH1A2 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA001566ED01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O94788 |
| Immunogen | Recombinant Human Retinal dehydrogenase 2 protein (1-518AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Retinal dehydrogenase 2 (RALDH 2) (RaLDH2) (EC 1.2.1.36) (Aldehyde dehydrogenase family 1 member A2) (Retinaldehyde-specific dehydrogenase type 2) (RALDH(II)), ALDH1A2, RALDH2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | ALDH1A2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |