



ARG1 Antibody, HRP conjugated

Product Code	CSB-PA002005DB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05089
Immunogen	Recombinant Human Arginase-1 protein (1-322AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Arginase-1 (EC 3.5.3.1) (Liver-type arginase) (Type I arginase), ARG1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.