



ARPC1B Antibody, HRP conjugated

Product Code	CSB-PA002126LB01HU
Abbreviation	Actin-related protein 2/3 complex subunit 1B
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15143
Immunogen	Recombinant Human Actin-related protein 2/3 complex subunit 1B protein (16-316AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Functions as component of the Arp2/3 complex which is involved in regulation of actin polymerization and together with an activating nucleation-promoting factor (NPF) mediates the formation of branched actin networks.
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin-related protein 2/3 complex subunit 1B (Arp2/3 complex 41 kDa subunit) (p41-ARC), ARPC1B, ARC41
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARPC1B