



# ARPC1B Antibody, FITC conjugated

Product Code	CSB-PA002126LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15143
Immunogen	Recombinant Human Actin-related protein 2/3 complex subunit 1B protein (16-316AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin-related protein 2/3 complex subunit 1B (Arp2/3 complex 41 kDa subunit) (p41-ARC), ARPC1B, ARC41
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARPC1B