



AVPR2 Antibody

Product Code	CSB-PA002470LA01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32307
Immunogen	Recombinant Pig Vasopressin V2 receptor protein (292-370AA)
Raised In	Rabbit
Species Reactivity	Pig
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vasopressin V2 receptor (V2R) (AVPR V2) (Antidiuretic hormone receptor) (Renal-type arginine vasopressin receptor), AVPR2
Immunogen Species	Sus scrofa (Pig)
Research Area	Others
Target Names	AVPR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.