




BET1 Antibody

Product Code	CSB-PA002670LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15155
Immunogen	Recombinant Human BET1 homolog protein (1-83AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BET1 homolog (hBET1) (Golgi vesicular membrane-trafficking protein p18), BET1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	BET1

Image



Western Blot
Positive WB detected in: K562 whole cell lysate,
Rat brain tissue
All lanes: BET1 antibody at 3.4µg/ml
Secondary
Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 14, 10 kDa
Observed band size: 14 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.