



C1QA Antibody, Biotin conjugated

Product Code	CSB-PA003637LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02745
Immunogen	Recombinant Human Complement C1q subcomponent subunit A protein (26-125AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement C1q subcomponent subunit A, C1QA
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	C1QA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.