




CCZ1 Antibody

Product Code	CSB-PA004156LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P86791
Immunogen	Recombinant Human Vacuolar fusion protein CCZ1 homolog protein (201-482AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar fusion protein CCZ1 homolog, CCZ1, C7orf28A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CCZ1

Image



Western Blot

Positive WB detected in: Rat liver tissue, Rat lung tissue

All lanes: CCZ1 antibody at 2.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 56 kDa

Observed band size: 56 kDa