



CCZ1 Antibody, HRP conjugated

Product Code	CSB-PA004156LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P86791
Immunogen	Recombinant Human Vacuolar fusion protein CCZ1 homolog protein (201-482AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar fusion protein CCZ1 homolog, CCZ1, C7orf28A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CCZ1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.