



CACNG3 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA004417LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O60359 |
| Immunogen | Recombinant Human Voltage-dependent calcium channel gamma-3 subunit protein (235-315AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Voltage-dependent calcium channel gamma-3 subunit (Neuronal voltage-gated calcium channel gamma-3 subunit) (Transmembrane AMPAR regulatory protein gamma-3) (TARP gamma-3), CACNG3 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | CACNG3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |