



CCR10 Antibody, FITC conjugated

Product Code	CSB-PA004840LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46092
Immunogen	Recombinant Human C-C chemokine receptor type 10 protein (314-362AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-C chemokine receptor type 10 (C-C CKR-10) (CC-CKR-10) (CCR-10) (G-protein coupled receptor 2), CCR10, GPR2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CCR10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.