




CD1C Antibody

Product Code	CSB-PA004891LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29017
Immunogen	Recombinant Human T-cell surface glycoprotein CD1c protein (18-302AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:500-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell surface glycoprotein CD1c (CD antigen CD1c), CD1C
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD1C

Image



Western Blot

Positive WB detected in: Mouse heart tissue


All lanes: CD1C antibody at 3μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 38 kDa

Observed band size: 38 kDa



Immunofluorescent analysis of A549 cells using CSB-PA004891LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)