



CD1C Antibody, FITC conjugated

Product Code	CSB-PA004891LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29017
Immunogen	Recombinant Human T-cell surface glycoprotein CD1c protein (18-302AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell surface glycoprotein CD1c (CD antigen CD1c), CD1C
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD1C