



CD3E Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA004931HB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P07766 |
| Immunogen | Recombinant Human T-cell surface glycoprotein CD3 epsilon chain protein (23-126AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | T-cell surface glycoprotein CD3 epsilon chain (T-cell surface antigen T3/Leu-4 epsilon chain) (CD antigen CD3e), CD3E, T3E |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | CD3E |