



CD3E Antibody, Biotin conjugated

Product Code	CSB-PA004931HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07766
Immunogen	Recombinant Human T-cell surface glycoprotein CD3 epsilon chain protein (23-126AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell surface glycoprotein CD3 epsilon chain (T-cell surface antigen T3/Leu-4 epsilon chain) (CD antigen CD3e), CD3E, T3E
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD3E
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.