





CEACAM6 Antibody

Product Code	CSB-PA005166LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40199
Immunogen	Recombinant Human Carcinoembryonic antigen-related cell adhesion molecule 6 protein (35-320AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carcinoembryonic antigen-related cell adhesion molecule 6 (Non-specific crossreacting antigen) (Normal cross-reacting antigen) (CD antigen CD66c), CEACAM6, NCA
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	CEACAM6

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA005166LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA005166LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.