



CH25H Antibody, Biotin conjugated

Product Code	CSB-PA005307LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95992
Immunogen	Recombinant Human Cholesterol 25-hydroxylase protein (142-272AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cholesterol 25-hydroxylase (EC 1.14.99.38) (Cholesterol 25-monooxygenase) (h25OH), CH25H
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CH25H