



chd Antibody, Biotin conjugated

Product Code	CSB-PA005377LD01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O57472
Immunogen	Recombinant Zebrafish Chordin protein (21-421AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chordin (Protein chordino), chd, chrd
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Others
Target Names	chd
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.