



# CLOCK Antibody, HRP conjugated

Product Code	CSB-PA005574EB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15516
Immunogen	Recombinant Human Circadian locomoter output cycles protein kaput protein (577-846AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Circadian locomoter output cycles protein kaput (hCLOCK) (EC 2.3.1.48) (Class E basic helix-loop-helix protein 8) (bHLHe8), CLOCK, BHLHE8 KIAA0334
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CLOCK