



CPB2 Antibody, Biotin conjugated

Product Code	CSB-PA005884LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96IY4
Immunogen	Recombinant Human Carboxypeptidase B2 protein (225-372AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carboxypeptidase B2 (EC 3.4.17.20) (Carboxypeptidase U) (CPU) (Plasma carboxypeptidase B) (pCPB) (Thrombin-activable fibrinolysis inhibitor) (TAFI), CPB2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CPB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.