



CR1 Antibody, HRP conjugated

Product Code	CSB-PA005932EB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17927
Immunogen	Recombinant Human Complement receptor type 1 protein (901-1095AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement receptor type 1 (C3b/C4b receptor) (CD antigen CD35), CR1, C3BR
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.