



Ctgf Antibody, Biotin conjugated

Product Code	CSB-PA006147LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29268
Immunogen	Recombinant Mouse Connective tissue growth factor protein (26-344AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Connective tissue growth factor (CCN family member 2) (Hypertrophic chondrocyte-specific protein 24) (Protein FISP-12), Ctgf, Ccn2 Fisp-12 Fisp12 Hcs24
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Ctgf
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.