



# CUTA Antibody, Biotin conjugated

Product Code	CSB-PA006226LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60888
Immunogen	Recombinant Human Protein CutA protein (33-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein CutA (Acetylcholinesterase-associated protein) (Brain acetylcholinesterase putative membrane anchor), CUTA, ACHAP C6orf82
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CUTA