



CUTA Antibody, Biotin conjugated

Product Code	CSB-PA006226LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60888
Immunogen	Recombinant Human Protein CutA protein (33-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein CutA (Acetylcholinesterase-associated protein) (Brain acetylcholinesterase putative membrane anchor), CUTA, ACHAP C6orf82
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CUTA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.