



DNASE1 Antibody, FITC conjugated

Product Code	CSB-PA007049LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24855
Immunogen	Recombinant Human Deoxyribonuclease-1 protein (78-231AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deoxyribonuclease-1 (EC 3.1.21.1) (Deoxyribonuclease I) (DNase I) (Dornase alfa), DNASE1, DNL1 DRNI
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DNASE1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.